

## NEAS Rabbit mAb

Catalog No: #52454

Package Size: #52454-1 50ul #52454-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NEAS Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-7C6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human NEAS
Conjugates	Unconjugated
Modification	Unmodification
Other Names	NEAS; EIEE5; SPTA2
Accession No.	Swiss-Prot:Q13813GeneID:6709
Uniprot	Q13813
GeneID	6709
Calculated MW	Calculated MW: 285 kDa; Observed MW: 285 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000;

## Images



Western blot detection of NEAS in mouse brain and rat brain lysates using NEAS antibody. Predicted band size: 285 kDa. Observed band size: 285 kDa.

## Background

---

Swiss-Prot Acc.Q13813.Fodrin, which seems to be involved in secretion, interacts with calmodulin in a calcium-dependent manner and is thus candidate for the calcium-dependent movement of the cytoskeleton at the membrane.

Regions	Feature key	Position(s)	Description
---------	-------------	-------------	-------------

---

Note: This product is for in vitro research use only