

NeuN Rabbit mAb

Catalog No: #52458



Package Size: #52458-1 50ul #52458-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NeuN Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S09-7H2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human NeuN
Conjugates	Unconjugated
Modification	Unmodification
Other Names	FOX3; NEUN; FOX-3; HRNBP3
Accession No.	Swiss-Prot:A6NFN3GeneID:
Uniprot	A6NFN3
Calculated MW	Calculated MW: 34 kDa; Observed MW: 46-55 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/1000;

Images



Western blot detection of NeuN in Mouse brain lysates using NeuN Rabbit mAb(1:1000 diluted).Predicted band size:34KDa.Observed band size: 46-55KDa.



Western blot detection of NeuN in Rat Brain lysates using NeuN Rabbit mAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:46-55kDa.

Background

Swiss-Prot Acc.A6NFN3.Pre-mRNA alternative splicing regulator. Regulates alternative splicing of RBFOX2 to enhance the production of mRNA species that are targeted for nonsense-mediated decay (NMD).

Note: This product is for in vitro research use only