

Nova1 Rabbit mAb

Catalog No: #52469

Package Size: #52469-1 50ul #52469-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Nova1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-8K6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Nova1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Nova-1
Accession No.	Swiss-Prot:P51513GenelD:
Uniprot	P51513
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/2000;

Images



Western blot detection of Nova1 in Rat lung cell lysates using Nova1 antibody. Observed band size:52kDa.

Background

Swiss-Prot Acc.P51513.May regulate RNA splicing or metabolism in a specific subset of developing neurons.

Note: This product is for in vitro research use only