

## ORP150 Rabbit mAb

Catalog No: #52480

Package Size: #52480-1 50ul #52480-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ORP150 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-4K7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ORP150
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Grp170; HSP12A; ORP150; GRP-170; ORP-150
Accession No.	Swiss-Prot:Q9Y4L1GeneID:
Uniprot	Q9Y4L1
Calculated MW	Calculated MW: 111 kDa; Observed MW: 150 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000-1/5000;

## Images



Western blot detection of ORP150 in 293T, MCF-7 using ORP150 Rabbit mAb. Predicted band size: 111kDa. Observed band size: 150kDa.

## Background

Swiss-Prot Acc.Q9Y4L1. Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular

chaperone and participate in protein folding.

---

Note: This product is for in vitro research use only