

p40 Rabbit mAb

Catalog No: #52482



Package Size: #52482-1 50ul #52482-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	p40 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-8K2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human p40
Conjugates	Unconjugated
Modification	Unmodification
Other Names	p40; RAB9p40; RABEPK; RGD1310612; RP11 65N13.1; RP23-446N16.2;
Accession No.	Swiss-Prot:Q7Z6M1GeneID:
Uniprot	Q7Z6M1
Calculated MW	Calculated MW: 41 kDa; Observed MW: 40 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of p40 in K562,HeLa cell lysates using p40 Rabbit mAb(1:1000 diluted).Predicted band size:41kDa.Observed band size:40kDa.

Background

Swiss-Prot Acc.Q7Z6M1.Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

Note: This product is for in vitro research use only