

PFKFB3 Rabbit mAb

Catalog No: #52501

Package Size: #52501-1 50ul #52501-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PFKFB3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-7A4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PFKFB3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	PFK2; IPFK2; iPFK-2
Accession No.	Swiss-Prot:Q16875GeneID:5209
Uniprot	Q16875
GeneID	5209
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000;

Images



Western blot detection of PFKFB3 in K562 cell lysates using PFKFB3 Rabbit mAb(1:1000 diluted).Predicted band size:60kDa.Observed band size:60kDa.

Background

Swiss-Prot Acc.Q16875.Synthesis and degradation of fructose 2,6-bisphosphate.

Note: This product is for in vitro research use only