

## PFKL Rabbit mAb

Catalog No: #52502

Package Size: #52502-1 50ul #52502-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PFKL Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-7A4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human PFKL
Conjugates	Unconjugated
Modification	Unmodification
Other Names	6-phosphofructokinase type B;PFK-B;Phosphohexokinase
Accession No.	Swiss-Prot:P17858GenelD:5211
Uniprot	P17858
GenelD	5211
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000;

## Images



Western blot detection of PFKL in C6,3T3 cell lysates using PFKL Rabbit mAb(1:1000 diluted).Predicted band size:85kDa.Observed band size:85kDa.

## Background

---

Swiss-Prot Acc.P17858. Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of glycolysis.

---

Note: This product is for in vitro research use only