

PPM1G Rabbit mAb

Catalog No: #52519

Package Size: #52519-1 50ul #52519-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PPM1G Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S05-8C4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PPM1G
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Protein phosphatase 1G;PP2C-gamma;Protein phosphatase magnesium-dependent 1 gamma
Accession No.	Swiss-Prot:O15355GeneID:5496
Uniprot	O15355
GeneID	5496
Calculated MW	Calculated MW: 59 kDa; Observed MW: 75 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000-1/10000;

Images



Western blot detection of PPM1G in K562,Rat Brain,C6,3T3,Hela cell lysates using PPM1G Rabbit mAb(1:1000 diluted).Predicted band size:59kDa.Observed band size:75kDa.

Background

Swiss-Prot Acc.O15355.Protein phosphatase 1G

Note: This product is for in vitro research use only