

PKA RI alpha Rabbit mAb

Catalog No: #52529



Package Size: #52529-1 50ul #52529-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PKA RI alpha Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-8A5
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Protein Kinase A regulatory subunit I alpha
Conjugates	Unconjugated
Modification	Unmodification
Other Names	CAR; CNC; CNC1; PKR1; TSE1; ADOHR; PPNAD1; PRKAR1; ACRDYS1
Accession No.	Swiss-Prot:P10644GenelD:5573
Uniprot	P10644
GenelD	5573
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000

Images



Western blot detection of Protein Kinase A regulatory subunit I alpha/PRKAR1A in C6,Rat Brain,Jurkat cell lysates using Protein Kinase A regulatory subunit I alpha/PRKAR1A Rabbit mAb(1:500 diluted).Predicted band size:43kDa.Observed band size:43kDa.

Background

Swiss-Prot Acc.P10644.Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells.

Note: This product is for in vitro research use only