

RSK1 p90 Rabbit mAb

Catalog No: #52563

Package Size: #52563-1 50ul #52563-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | RSK1 p90 Rabbit mAb |
| Host Species | Recombinant Rabbit |
| Clonality | Monoclonal antibody |
| Clone No. | S02-2D2 |
| Isotype | Rabbit IgG |
| Purification | Affinity Purified |
| Applications | WB IHC |
| Species Reactivity | Human,Mouse,Rat |
| Immunogen Description | A synthetic peptide of human RSK1 p90 |
| Conjugates | Unconjugated |
| Modification | Unmodification |
| Other Names | RSK; HU-1; RSK1; p90Rsk; MAPKAPK1A |
| Accession No. | Swiss-Prot:Q15418GeneID:6195 |
| Uniprot | Q15418 |
| GeneID | 6195 |
| Calculated MW | Calculated MW: 83 kDa; Observed MW: 90 kDa |
| Concentration | 0.3 mg/ml |
| Formulation | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB: 1/1000; IHC: 1/20;

Images



Western blot detection of RSK1 p90 in HeLa cell lysates using RSK1 p90 Rabbit mAb(1:500 diluted).Predicted band size:83kDa.Observed band size:90kDa.



Immunohistochemistry of RSK1 p90 in paraffin-embedded Human Cholangiocarcinoma using RSK1 p90 Rabbit mAb at dilution 1/50

Background

Swiss-Prot Acc.Q15418.Serine/threonine-protein kinase that acts downstream of ERK (MAPK1/ERK2 and MAPK3/ERK1) signaling and mediates mitogenic and stress-induced activation of the transcription factors CREB1, ETV1/ER81 and NR4A1/NUR77, regulates translation through RPS6 and EIF4B phosphorylation, and mediates cellular proliferation, survival, and differentiation by modulating mTOR signaling and repressing pro-apoptotic function of BAD and DAPK1. In fibroblast, is required for EGF-stimulated phosphorylation of CREB1, which results in the subsequent transcriptional activation of several immediate-early genes. In response to mitogenic stimulation (EGF and PMA), phosphorylates and activates NR4A1/NUR77 and ETV1/ER81 transcription factors and the cofactor CREBBP. Upon insulin-derived signal, acts indirectly on the transcription regulation of several genes by phosphorylating GSK3B at Ser-9; and inhibiting its activity. Phosphorylates RPS6 in response to serum or EGF via an mTOR-independent mechanism and promotes translation initiation by facilitating assembly of the pre-initiation complex. In response to insulin, phosphorylates EIF4B, enhancing EIF4B affinity for the EIF3 complex and stimulating cap-dependent translation. Is involved in the mTOR nutrient-sensing pathway by directly phosphorylating TSC2 at Ser-1798;, which potently inhibits TSC2 ability to suppress mTOR signaling, and mediates phosphorylation of RPTOR, which regulates mTORC1 activity and may promote rapamycin-sensitive signaling independently of the PI3K/AKT pathway. Mediates cell survival by phosphorylating the pro-apoptotic proteins BAD and DAPK1 and suppressing their pro-apoptotic function. Promotes the survival of hepatic stellate cells by phosphorylating CEBPB in response to the hepatotoxin carbon tetrachloride (CCl4). Mediates induction of hepatocyte proliferation by TGFA through phosphorylation of CEBPB. Is involved in cell cycle regulation by phosphorylating the CDK inhibitor CDKN1B, which promotes CDKN1B association with 14-3-3 proteins and prevents its translocation to the nucleus and inhibition of G1 progression. Phosphorylates EPHA2 at Ser-897;, the RPS6KA-EPHA2 signaling pathway controls cell migration (PubMed:26158630).

Note: This product is for in vitro research use only