

Sall4 Rabbit mAb

Catalog No: #52565

Package Size: #52565-1 50ul #52565-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Sall4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-6H7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Sall4
Conjugates	Unconjugated
Modification	Unmodification
Other Names	DRRS; HSAL4; ZNF797
Accession No.	Swiss-Prot:Q9UJQ4GeneID:57167
Uniprot	Q9UJQ4
GeneID	57167
Calculated MW	Calculated MW: 112 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of Sall4 in Daudi using Sall4 Rabbit mAb. Predicted band size: 112kDa. Observed band size: 112kDa.

Background

Swiss-Prot Acc.Q9UJQ4. Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.

Note: This product is for in vitro research use only