

SLCO5A1 Rabbit mAb

Catalog No: #52582

Package Size: #52582-1 50ul #52582-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SLCO5A1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-5J4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human SLCO5A1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	OATPJ; OATP-J; OATP5A1; OATPRP4; OATP-RP4; SLC21A15
Accession No.	Swiss-Prot:Q9H2Y9GeneID:81796
Uniprot	Q9H2Y9
GeneID	81796
Calculated MW	Calculated MW: 92 kDa; Observed MW: 92 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of SLCO5A1 in K562, HeLa cell lysates using SLCO5A1 Rabbit mAb (1:1000 diluted). Predicted band size: 92kDa. Observed band size: 92kDa.

Background

Swiss-Prot Acc.Q9H2Y9.Solute carrier organic anion transporter family member 5A1

Note: This product is for in vitro research use only