

Somatostatin Receptor 3 Rabbit mAb

Catalog No: #52587

Package Size: #52587-1 50ul #52587-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Somatostatin Receptor 3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-2G8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Somatostatin Receptor 3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	SS3R; SS3-R; SS-3-R; SSR-28
Accession No.	Swiss-Prot:P32745GenelD:6753
Uniprot	P32745
GenelD	6753
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000-1/10000;

Images



Western blot detection of Somatostatin Receptor 3/SSTR3 in Jurkat cell lysates using Somatostatin Receptor 3/SSTR3 Rabbit mAb(1:500 diluted). Predicted band size:46kDa.Observed band size:46kDa.

Background

Swiss-Prot Acc.P32745.Receptor for somatostatin-14 and -28. This receptor is coupled via pertussis toxin sensitive G proteins to inhibition of adenylyl cyclase.

Note: This product is for in vitro research use only