

Tensin 4 Rabbit mAb

Catalog No: #52624

Package Size: #52624-1 50ul #52624-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Tensin 4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-1F3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human TNS4
Conjugates	Unconjugated
Modification	Unmodification
Other Names	CTEN; C-terminal tensin-like protein
Accession No.	Swiss-Prot:Q8IZW8GeneID:84951
Uniprot	Q8IZW8
GeneID	84951
Calculated MW	Calculated MW: 77 kDa; Observed MW: 90 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of TNS4/CTEN in K562 cell lysates using TNS4/CTEN Rabbit mAb(1:1000 diluted).Predicted band size:77kDa.Observed band size:90kDa.

Background

Swiss-Prot Acc.Q8IZW8.May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton . May promote apoptosis, via its cleavage by caspase-3.

Note: This product is for in vitro research use only