

## NUP62CL Rabbit mAb

Catalog No: #52682

Package Size: #52682-1 50ul #52682-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NUP62CL Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-1A5
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human NUP62CL
Conjugates	Unconjugated
Modification	Unmodification
Other Names	NUP62L; Nucleoporin-62 C-terminal-like protein
Accession No.	Swiss-Prot:Q9H1M0GeneID:54830
Uniprot	Q9H1M0
GeneID	54830
Calculated MW	Calculated MW:21 kDa,Observed MW:21 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of NUP62CL in HeLa,Ramos cell lysates using NUP62CL Rabbit mAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.

## Background

---

This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore complexes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]

---

Note: This product is for in vitro research use only