

Bcl2 L14 Rabbit mAb

Catalog No: #52721

Package Size: #52721-1 50ul #52721-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Bcl2 L14 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-3D8
Isotype	IgG
Purification	Affinity Purified
Applications	WB IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Bcl G/BCL2L14
Conjugates	Unconjugated
Modification	Unmodification
Other Names	BCLG
Accession No.	Swiss-Prot:Q9BZR8GeneID:79370
Uniprot	Q9BZR8
GeneID	79370
Calculated MW	Calculated MW:37 kDa,Observed MW:43 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


WB: 1/1000 ICC/IF: 1/50

Images



Western blot detection of Bcl G/BCL2L14 in HeLa,A549,HL-60,U251,U87-MG cell lysates using Bcl G/BCL2L14 Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:43kDa.

Immunocytochemistry of Bcl G/BCL2L14 (green) in hela using Bcl G/BCL2L14 Rabbit mAb at dilution 1/50, and DAPI(blue)



Background

Plays a role in apoptosis.

Note: This product is for in vitro research use only