

## Bid Rabbit mAb

Catalog No: #52725

Package Size: #52725-1 50ul #52725-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Bid Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-2F0
Isotype	IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Bid
Conjugates	Unconjugated
Modification	Unmodification
Other Names	FP497; MGC15319; MGC42355; BID
Accession No.	Swiss-Prot:P55957GenelD:637
Uniprot	P55957
GenelD	637
Calculated MW	Calculated MW:22 kDa,Observed MW:22 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/1000; IHC: 1/50;

## Images



Western blot detection of Bid in HeLa,A549,HL-60 using Bid Rabbit mAb(1:1000 diluted)



Immunohistochemistry of Bid in paraffin-embedded liver cancer tissue using Bid Rabbit mAb at dilution 1/50

## Background

Induces caspases and apoptosis (PubMed:14583606).

Counters the protective effect of BCL2 (By similarity).

Note: This product is for in vitro research use only