

DLX3 Rabbit mAb

Catalog No: #52750

Package Size: #52750-1 50ul #52750-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DLX3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-0C2
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human DLX3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AI4; TDO
Accession No.	Swiss-Prot:O60479GeneID:1747
Uniprot	O60479
GeneID	1747
Calculated MW	Calculated MW:32 kDa,Observed MW:32 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of DLX3 in A549,U2OS cell lysates using DLX3 Rabbit mAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32kDa.

Background

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Note: This product is for in vitro research use only