

## DNAJC12 Rabbit mAb

Catalog No: #52751

Package Size: #52751-1 50ul #52751-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DNAJC12 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-3E8
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human DNAJC12
Conjugates	Unconjugated
Modification	Unmodification
Other Names	JDP1; HPANBH4
Accession No.	Swiss-Prot:Q9UKB3GeneID:56521
Uniprot	Q9UKB3
GeneID	56521
Calculated MW	Calculated MW:23 kDa,Observed MW:23 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of DNAJC12 in HeLa,A549,HL-60 cell lysates using DNAJC12 Rabbit mAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.

## Background

---

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

---

Note: This product is for in vitro research use only