

## HDLBP Rabbit mAb

Catalog No: #52770

Package Size: #52770-1 50ul #52770-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HDLBP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-3H2
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human HDLBP
Conjugates	Unconjugated
Modification	Unmodification
Other Names	HBP; VGL; PRO2900
Accession No.	Swiss-Prot:Q00341GeneID:3069
Uniprot	Q00341
GeneID	3069
Calculated MW	Calculated MW:142 kDa,Observed MW:142 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/2000-1/10000

## Images



Western blot detection of HDLBP in HeLa,A549,HL-60,U251,U87-MG cell lysates using HDLBP Rabbit mAb(1:1000 diluted).Predicted band size:142kDa.Observed band size:142kDa.

## Background

---

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

---

Note: This product is for in vitro research use only