

ITCH Rabbit mAb

Catalog No: #52785

Package Size: #52785-1 50ul #52785-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ITCH Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-416
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ITCH/AIP4
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AIF4; AIP4; ADMFD; NAPP1
Accession No.	Swiss-Prot:Q96J02GeneID:83737
Uniprot	Q96J02
GeneID	83737
Calculated MW	Calculated MW:103 kDa,Observed MW:103 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of ITCH/AIP4 in A549,HL-60 cell lysates using ITCH/AIP4 Rabbit mAb(1:1000 diluted).Predicted band size:103kDa.Observed band size:103kDa.

Background

Acts as an E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates?

Note: This product is for in vitro research use only