

Nectin 2 Rabbit mAb

Catalog No: #52817

Package Size: #52817-1 50ul #52817-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Nectin 2 Rabbit mAb
Clone No.	S06-8F3
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB,ICC/IF,IP,FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Nectin 2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	HVEB; PRR2; CD112; PVRL2; PVRR2
Accession No.	Swiss-Prot:Q92692GeneID:5819
Uniprot	Q92692
GeneID	5819
Calculated MW	62-75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:100 FC 1:50

Images



Western blot analysis of Nectin 2 expression in SW480 cell lysate.

Background

Probable cell adhesion protein.

Note: This product is for in vitro research use only