

Pleckstrin Rabbit mAb

Catalog No: #52840

Package Size: #52840-1 50ul #52840-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Pleckstrin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-7C5
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Pleckstrin
Conjugates	Unconjugated
Modification	Unmodification
Other Names	P47;Platelet 47 kDa protein
Accession No.	Swiss-Prot:P08567GenelD:5341
Uniprot	P08567
GenelD	5341
Calculated MW	Calculated MW:40 kDa,Observed MW:40 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of Pleckstrin in HeLa,HL-60 cell lysates using Pleckstrin Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:40kDa.

Background

Major protein kinase C substrate of platelets.

Note: This product is for in vitro research use only