

STAT6 Rabbit mAb

Catalog No: #52868

Package Size: #52868-1 50ul #52868-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	STAT6 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-0B5
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human STAT6
Conjugates	Unconjugated
Modification	Unmodification
Other Names	STAT6B; STAT6C; D12S1644; IL-4-STAT
Accession No.	Swiss-Prot:P42226GenelD:6778
Uniprot	P42226
GenelD	6778
Calculated MW	Calculated MW:94 kDa,Observed MW:110 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000

Images



Western blot detection of STAT6 in HeLa,U2OS cell lysates using STAT6 Rabbit mAb(1:1000 diluted).Predicted band size:94kDa.Observed band size:94kDa.

Background

Carries out a dual function: signal transduction and activation of transcription. Involved in IL4/interleukin-4- and IL3/interleukin-3-mediated signaling.

Note: This product is for in vitro research use only