

USF1 Rabbit mAb

Catalog No: #52885



Package Size: #52885-1 50ul #52885-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	USF1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-5E6
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human USF1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	UEF; FCHL; MLTF; FCHL1; MLTF1; HYPLIP1; bHLHb11
Accession No.	Swiss-Prot:P22415GenelD:7391
Uniprot	P22415
GenelD	7391
Calculated MW	Calculated MW:34 kDa,Observed MW:34 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000-1/10000

Images



Western blot detection of USF1 in HeLa,mouse Brain,Ramos using USF1 Rabbit mAb(1:1000 diluted)

Background

Transcription factor that binds to a symmetrical DNA sequence (E-boxes) (5'-CACGTG-3') that is found in a variety of viral and cellular promoters.

Note: This product is for in vitro research use only