

DAP kinase 3 Rabbit mAb

Catalog No: #52889

Package Size: #52889-1 50ul #52889-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DAP kinase 3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-7F5
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZIP Kinase
Conjugates	Unconjugated
Modification	Unmodification
Other Names	DLK; ZIP; ZIPK
Accession No.	Swiss-Prot:O43293GeneID:1613
Uniprot	O43293
GeneID	1613
Calculated MW	Calculated MW:53 kDa,Observed MW:53 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of ZIP Kinase in HeLa,A549,MCF-7,U251 cell lysates using ZIP Kinase Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.

Background

Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis. [provided by RefSeq, Jul 2008]

Note: This product is for in vitro research use only