

ZNF217 Rabbit mAb

Catalog No: #52890

Package Size: #52890-1 50ul #52890-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF217 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-4B1
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZNF217
Conjugates	Unconjugated
Modification	Unmodification
Other Names	ZABC1
Accession No.	Swiss-Prot:O75362GeneID:7764
Uniprot	O75362
GeneID	7764
Calculated MW	Calculated MW:115 kDa,Observed MW:115 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000-1/10000

Images



Western blot detection of ZNF217 in A549,HL-60,U251 cell lysates using ZNF217 Rabbit mAb(1:1000 diluted).Predicted band size:115kDa.Observed band size:115kDa.

Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Note: This product is for in vitro research use only