

IL-17D Polyclonal Antibody

Catalog No: #29667

Package Size: #29667-1 50ul #29667-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-17D Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IHC-p,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous IL-17D
Immunogen Description	Synthetic peptide from human protein at AA range: 141-190
Accession No.	Swiss-Prot#: Q8TAD2
Uniprot	Q8TAD2
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100

Note: This product is for in vitro research use only