

Rat N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit



Catalog No: #EK20248

Orders: order@signalwayantibody.com

Package Size: #EK20248-1 48T #EK20248-2 96T

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------|
| Product Name | Rat N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rt |

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.128ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only