

Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit

Catalog No: #EK20254

Orders: order@signalwayantibody.com

Package Size: #EK20254-1 48T #EK20254-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Accession No.	Swiss Prot:Q9CPV9
Uniprot	Q9CPV9

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.061ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only