

# Mouse Interleukin 4 Induced Protein 1 (IL4I1) ELISA Kit



Catalog No: #EK20270

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20270-1 48T #EK20270-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Interleukin 4 Induced Protein 1 (IL4I1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:24.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only