

# Rat Endothelial microparticle (EMP) ELISA Kit

Catalog No: #EK20273



Package Size: #EK20273-1 48T #EK20273-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Endothelial microparticle (EMP) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:11.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only