

Human Extracellular Matrix Protein 1 (ECM1) ELISA Kit



Catalog No: #EK20315

Orders: order@signalwayantibody.com

Package Size: #EK20315-1 48T #EK20315-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Extracellular Matrix Protein 1 (ECM1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range: 31.2-2000pg/mL

Sensitivity: 13.2pg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only