

# Rat Interferon Regulatory Factor 7 (IRF7) ELISA Kit

Catalog No: #EK20318



Package Size: #EK20318-1 48T #EK20318-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Rat Interferon Regulatory Factor 7 (IRF7) ELISA Kit |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Rat   |

## Application Details

|  |
|--|
| Detect Range:62.5-4000pg/mL  |
| Sensitivity:22.3pg/mL  |
| Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume: 50 &mu;L  |
| Assay Time:1-4.5h  |
| Detection wavelength:450 nm  |

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only