

# ZBTB2 Polyclonal Antibody

Catalog No: #44015



Package Size: #44015-1 50ul #44015-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZBTB2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB ELISA
Species Reactivity	Human,Rat,Mouse
Specificity	ZBTB2 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	ZBTB2 KIAA1483 ZNF437
Accession No.	SwissProt:Q8N680
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only