

CD63 Polyclonal Conjugated Antibody

Catalog No: #C44012



Package Size: #C44012-AF350 100ul #C44012-AF405 100ul #C44012-AF488 100ul
 #C44012-AF555 100ul #C44012-AF594 100ul #C44012-AF647 100ul
 #C44012-AF680 100ul #C44012-AF750 100ul #C44012-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CD63 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Human,Rat,Mouse
Specificity	CD63 Polyclonal Antibody detects endogenous levels of CD63 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from the Internal region of human CD63.
Accession No.	SwissProt:P08962
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only