

## Recombinant Human Homeobox protein SIX2(SIX2)

Catalog No: #AP89681

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Homeobox protein SIX2(SIX2)
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	Homeobox protein SIX2; OTTHUMP00000201649; Sine oculis homeobox (Drosophila) homolog 2; Sine oculis homeobox homolog 2 (Drosophila); Sine oculis homeobox homolog 2; SIX homeobox 2; Six2; SIX2_HUMAN
Accession No.	Q9NPC8
Uniprot	Q9NPC8
Target Species	Hu
Tag Info	N-terminal His-tagged
Target Sequence	1-291 full length protein MSMLPTFGFT QEQVACVCEV LQQGGNIERL GRFLWSPAC EHLHKNESVL KAKAVVAFHR GNFRELYKIL ESHQFSPHNH AKLQQLWLKA HYIEAEKLRG RPLGAVGKYR VRRKFPLPRS IWDGEETSYC FKEKRSRSLR EWYAHNPYPY PREKRELAEA TGLTTTQVSN WFKNRRQRDR AAEAKERENN ENSNSNSHNP LNGSGKSVLG SSEDEKTPSG TPDHSSSSPA LLLSPPPPGL PSLHSLGHPP GPSAVPVPVP GGGGADPLQH HHGLQDSILN PMSANLVDLG S
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. AP89681-2

## Application Details

&gt;85% (SDS-PAGE)

Note: This product is for in vitro research use only