

Recombinant Human Serine/threonine-protein kinase LATS1(LATS1) ,partial



Catalog No: #AP89694

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Serine/threonine-protein kinase LATS1(LATS1) ,partial
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	h-warts; Large tumor suppressor homolog 1; LATS large tumor suppressor homolog 1; LATS1; LATS1_HUMAN; Serine threonine protein kinase LATS1; Serine/threonine-protein kinase LATS1; WARTS; WARTS protein kinase; wts
Accession No.	O95835
Uniprot	O95835
Target Species	Hu
Tag Info	N-terminal His-tagged
Target Sequence	Partial
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> <p>AP89694-2</p>

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only