

## PSMD12 Polyclonal Antibody

Catalog No: #44017

Package Size: #44017-1 50ul #44017-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PSMD12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	PSMD12 Polyclonal Antibody detects endogenous levels of PSMD12 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human PSMD12.
Accession No.	Swiss-Prot#: O00232
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

## Images



Western Blot analysis of HepG2 cells using PSMD12 Polyclonal Antibody

Note: This product is for in vitro research use only