

# Mouse Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK20330



Package Size: #EK20330-1 48T #EK20330-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Immunoglobulin G (IgG) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range: 1.56-100 µg/mL

Sensitivity: 0.44 µg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only