

Mouse Osteomodulin (OMD) ELISA Kit

Catalog No: #EK20341



Package Size: #EK20341-1 48T #EK20341-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Osteomodulin (OMD) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000 pg/mL

Sensitivity:5.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only