

Rat Four And A Half LIM Domains Protein 1 (FHL1) ELISA Kit

Catalog No: #EK20334

Orders: order@signalwayantibody.com

Package Size: #EK20334-1 48T #EK20334-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat Four And A Half LIM Domains Protein 1 (FHL1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:25.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only