

YTHD2 Polyclonal Antibody

Catalog No: #44018



Package Size: #44018-1 50ul #44018-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	YTHD2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	YTHD2 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	YTHDF2 HGRG8
Accession No.	swiss-prot:Q9Y5A9
Calculated MW	63 kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

ELISA 1:5000-20000

Note: This product is for in vitro research use only