

SIA4A rabbit pAb

Catalog No: #44019



Package Size: #44019-1 50ul #44019-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SIA4A rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human, Mouse
Specificity	This antibody detects endogenous levels of SIA4A at Human/Mouse
Immunogen Description	Synthesized peptide derived from human SIA4A
Accession No.	Swiss-Prot#: Q11201
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1 μ g 500-2000

Note: This product is for in vitro research use only