

Recombinant Human ATP-binding cassette sub-family G member 2(ABCG2)

Catalog No: #AP89640

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89640-1 20ug #AP89640-2 100ug #AP89640-3 1mg

Description

Product Name	Recombinant Human ATP-binding cassette sub-family G member 2(ABCG2)
Host Species	Yeast
Immunogen Description	Expression Region: 557-630aa Sequence Info:Partial
Other Names	ABCG2; ABCP; BCRP; BCRP1; MXR; Broad substrate specificity ATP-binding cassette transporter ABCG2; ATP-binding cassette sub-family G member 2; Breast cancer resistance protein; CDw338; Mitoxantrone resistance-associated protein; Placenta-specific ATP-binding cassette transporter; Urate exporter; CD antigen CD338
Accession No.	Q9UNQ0
Target Species	Homo sapiens (Human)
SDS-PAGE MW	12.1 kDa
Tag Info	C-terminal 6xHis-Myc-tagged
Target Sequence	NLTTIASWLSWLQYFSIPRYGFTALQHNEFLGQNFPCPLNATGNNPCNYATCTGEEYLVKQGIDLSPWGLWK NH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only