

# Mouse Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit

Catalog No: #EK20336



Package Size: #EK20336-1 48T #EK20336-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit
Brief Description	ELISA Kit
Species Reactivity	Ms

## Application Details

Detect Range:3.12-200 mIU/mL

Sensitivity:1.42 mIU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only