

# Human CMP-N-acetylneuraminate-beta-galactosamide-alpha- 2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit

Catalog No: #EK20337

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20337-1 48T #EK20337-2 96T

## Description

Product Name	Human CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit
Brief Description	ELISA Kit
Species Reactivity	Human

## Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only