

ASAH3L Polyclonal Antibody

Catalog No: #44025



Package Size: #44025-1 50ul #44025-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ASAH3L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB;IHC-p;IF;ELISA
Species Reactivity	Human,Mouse,Monkey
Specificity	ASAH3L Polyclonal Antibody detects endogenous levels of ASAH3L protein.
Immunogen Description	Synthesized peptide derived from ASAH3L .
Accession No.	Swiss-Prot#: Q5QJU3
Calculated MW	33kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000

IHC-p: 1:100-1:300

Note: This product is for in vitro research use only