

SGPP1 rabbit pAb

Catalog No: #44044



Package Size: #44044-1 50ul #44044-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SGPP1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Synthesized peptide derived from human SGPP1
Accession No.	Swiss-Prot#: Q9BX95
Concentration	1mg/ml
Storage	Store at -20°C/1 year

Application Details

WB 1:1000 500-2000

Note: This product is for in vitro research use only