

# Anti-PD-1 [Nivolumab]

Catalog No: #52899



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Anti-PD-1 [Nivolumab]
Clonality	Monoclonal
Isotype	IgG4
Applications	Functional Assay
Specificity	Humanized Monoclonal Antibody
Other Names	Programmed cell death protein 1; Pdcd1; CD279
Formulation	0.01M PBS, pH 7.2
Storage	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at $\beta$ 80°C. Avoid repeated freeze/thaw cycles.

## Background

Nivolumab biosimilar is a human monoclonal antibody that blocks the interaction between PD-1 and its ligands, PD-L1 and PD-L2.

Note: This product is for in vitro research use only